



Self-Check-In Software for BYO Tablet on Stand



Image supplied by <https://www.bosstab.com/au/>

BOSSTAB

Why you need AutoMed's Self-Check-In Software for BYO Tablet

- ✓ AutoMed has released an alternative simple Self-Check-In Service for our Complete Services Solutions Packaged Clinics.
- ✓ The software is ideal for Clinics that already have Check-In hardware (i.e. a Tablet on a Stand) and would like a simple Check-in Service in lieu of the AutoMed Kios
- ✓ The software easily can be installed onto the Tablet.
- ✓ The software is available to Individual Clinics and also to Clinics Multiple Locations on a Shared Server.



What does the Clinic Need to provide?

- ✓ Clinics will need to provide their own Hardware - Tablet Device/s + Secure stand/s.
- ✓ Clinics without existing Hardware may purchase a Tablet + Secure Floor Stand.

i.e. Bosstab Elite Evo Floor Stand with Lockable Tablet Holder

<https://www.bosstab.com/au/tablet-stands/elite/>

Compatible with:

- 10.2" iPad Gen 7/8/9 | 9.7" iPad Air Gen 2 / iPad Pro | 12.9" iPad Pro Gen 3/4/5/6.
- 10.4" Samsung Galaxy Tab A7 | 10.5" Samsung Galaxy Tab A8 | 11" Samsung Galaxy Tab S7/S8 | 12.4" Samsung Galaxy Tab S7+/S8+ .
- 13" Microsoft Surface Pro Gen 8/9.

Settings

- ✓ Clinics will need to install “Kiosk-Mode” software on their tablet to make the device available as a kiosk.
- ✓ This app needs to restrict user from accidentally or purposely navigating to the home-screen.
 - ❖ iOS kiosk mode is a restrictive mode that locks down an iPad to a single app or a specific set of apps to run in the foreground.
 - i.e. Kiosk Pro Lite, Hexnode etc.
 - ❖ Android Kiosk Mode is a lock down mechanism that enables IT admins to run a single app or a pre-defined set of apps on Android smartphones and tablets, allowing these devices to function as kiosks.
 - i.e. Fully Kiosk, Android Kiosk, Hexnode, Boss tab etc.

Setup

- ✓ The Self-Check-In software will be in the form of a dedicated URL that will be supplied by AutoMed.
- ✓ The URL will be supplied to each clinic or location along with a once off activation PIN to install the software on their own hardware (Tablet).
- ✓ For Security purposes, the Clinic will need to advise AutoMed the Tablet IP address of the Tablet Device/s.

Cost

- ✓ The software is available at a cost per clinic or per location.
 - i.e. Packaged Clinics will change to a Complete Solutions Services Package with Self-Check-In
 - Medium Clinic \$415/mth.
 - Multiple locations on a shared server - Initial Package increase to packaged cost (as above) + \$50/mth per additional location.
- ✓ There are no Set-up or Configuration Fees.
- ✓ BYO Device + Stand is required.
- ❖ All pricing is advised exclusive of GST



Bosstab Elite Evo Floor Stand with Lockable Tablet Holder <https://www.bosstab.com/au/tablet-stands/elite/>
Compatible with o 10.2" iPad Gen 7/8/9 | 9.7" iPad Air Gen 2 / iPad Pro | 12.9" iPad Pro Gen 3/4/5/6
10.4" Samsung Galaxy Tab A7 | 10.5" Samsung Galaxy Tab A8 | 11" Samsung Galaxy Tab S7/S8
12.4" Samsung Galaxy Tab S7+/S8+ | 13" Microsoft Surface Pro Gen 8/9.


Patient Journey

- ✓ The following slides are the screens that the patient will see during the Self-Check-In process.

AutoMed Self-Check-In




Tap to Begin



Do you have an appointment today?

AutoMed Demo Clinic
(AutoMed Test Clinic)


Start Again



Do you have an appointment today?

Have you visited “Clinic Name” before?

AutoMed Demo Clinic
(AutoMed Test Clinic) Start Again



Have you visited AutoMed Test Clinic before?

No Yes

Enter Patient Surname

AutoMed Demo Clinic
(AutoMed Test Clinic) Start Again

Patient Surname

Patient Surname

Next

Enter Date of Birth:

AutoMed Demo Clinic
(AutoMed Test Clinic) Start Again


Date of Birth

DD MM YYYY

Back Next

Is this your appointment?


AutoMed Demo Clinic
(AutoMed Test Clinic) Start Again



Is this your appointment?
Standard appt.
at 03:20 PM with Dr James Frederick

Is this your Address?


AutoMed Demo Clinic
(AutoMed Test Clinic) Start Again



Is this your address?
114c Watton Street, Werribee, 3030

You have checked in - estimated wait time advised


AutoMed Demo Clinic
(AutoMed Test Clinic) Start Again




You have checked in
Estimated waiting time 20 minutes

Close

If the patient does not have an appointment or has not visited the clinic before the screen will advise:




If at any time patient details are unable to be located, a pop-up will advise:



Troubleshooting

- ✓ The following screen advises that the IP address has not been whitelisted.
- ✓ Please email immediately with the Tablet IP address to be loaded against the clinics profile.
- ✓ alison@automedsystems.com.au or support@automedsystems.com.au



Further Information

- ✓ Please contact the Sales Team directly to request the software as part of your package.
- ✓ Please direct any questions to sales@automedsystems.com.au or alison@automedsystems.com.au
- ✓ Subject - Clinic Name - Self-Check-In software for BYO Tablet